World First Slim Cassette Type Digital Mammo. Upgrade Solution





Contents

- **01** RoseM
 - What is RoseM?
 - Site Installation Photos
 - How Does RoseM Operate?
 - Going Digital wit RoseM
- **02** RoseM's Benefits
- **03** Performance
 - Image Quality
 - Image Comparison

- **04** Core Technologies
 - TRUVIEWTM ART
 - Lossless AED
 - RConsole1 X-ray Imaging Software
- **05** Market Performance
 - Proven Market Performance in Korea
 - RoseM Marketing in Korean Market
- **06** Further Info
 - Testimonials
 - Q&A
- **07** Specifications
 - Specifications
 - AEC Compatibility



Company Information

Company Name

DRTECH Corp.

Establish Date

Apr. 11, 2000

HQ

R&D HQ - Korea Factory - Korea Employees


152 (Nov. 2016)




General Radiography
Mammography
Veterinary
EOD/NDT
Mobile System
Portable Accessories


Global Operations

DRTECH Shanghai Co., Ltd DRTECH America Service Center DRTECH R&D Office of America DRTECH Europe GmbH DRTECH Sales Office of Turkey











Product Portfolio

Exprimer

• EVS 3643

Wireless

• EVS 4343

Wired

• EVS 2430

Wireless/Wired

FLAATZ

- F600
- F560
- F760
- F500
- F750E
- F330



- RSM 1824C
- RSM 2430C
- RSM 1824TD
- RSM 2430TD
- RSM 2430UD Under Development

ECONSOLE1

General Radiography

• EConsole1 + Truview® ART

ECONSOLE1 VET

Veterinary

EConsole1 Vet + Truview[®] ART

Veterinary System

• CUVIC 56



Mammography

• RConsole1 + Truview® ART

Portable Solution

• PCP

Portable Console PC

•WPCS

Wireless Power Charging System



What is RoseM?

DRTECH's latest innovation that introduced the concept of 'Retrofit' to the mammography market! With RoseM, an Analog System is instantaneously upgraded to a digital system.



Easy Upgrade

Simply insert RoseM into the existing bucky of an analog system for an Easy Upgrade to digital



Excellent Image Quality

76um pixel size with **Truview**[®] **ART** image processing algorithm; ensures acquisition and processing of **high quality image quality** paralleling top-tier Digital Mammo. Systems



Slim Cassette Size

Perfect compatibility with analog mammo. systems the size of RoseM is equal to any Film or CR IP(Imaging Plate) cassettes available on the market.



Perfect AEC Calibration

Designed to be compatible with AEC sensor of all existing analog systems.





What is RoseM?

Proven and **Reliable** Technology Approved by Regulatory Authorities

RSM1824C










RSM2430C












What is RoseM?



Perfect AEC Compatibility with any Analog System

AEC Compatibility technology Patent registered and approval pending





RoseM Site Installation Photos









How does RoseM operate?

	Rose M	Other methods
DR Upgrade Method summary	Same as CR or Film cassette	Entire bucky change and x-ray signal integration for PC SW required
Changing part	None, only AEC Calibration needed like CR or Film cassette	Whole Bucky with DR Mammo detector or modify the bucky for cable connecting Electronic board for X-ray signal integration with detector Cables for interconnecting X-ray generator to Detector Cables for interconnecting X-ray generator to PC
Certification	CE, MFDS, FDA cleared	Not cleared, Needed to be checked
Compatibility	All analogue system	Only system maker can modify
Time for install	Same as CR install (Several hours)	Long time to modify all X-ray system



Going Digital with RoseM

01 Digital Upgrade Solution

"Reengineer your old Mammo system with digital benefits of RoseM"



02 New Digital Equipment

"Get a new Digital system with benefits of RoseM"



03 Mobile Health Buses

"High Resistance to vibration & External Force"





RoseM's Benefits

DOUBLE BENEFITS





Comments by leading radiologists & Breast Cancer Physicians on RoseM Image Quality

"The images are even better than Hologic."

"It deserves premium pricing as the image quality is excellent."

"You should promote its image samples to win in market."



RoseM's Benefits

Maintenance & Consumable cost savings

- · Elimination of consumables Film / Chemicals / Image Plates etc
- · Reduced maintenance and service costs

Fast patient outflow

(=increased no. of patients)

- · Reduced patient waiting time
- · Increased economical benefits for hospitals with increased number of patients

Increased Efficiency with Digitalization

- · Reduced of Film storage costs and space
- · Less number of Film maintenance personnels

RoseM



Consumables None

Maintenance Cost



RoseM





Screen



PACS





Only

10 sec!

Diagnosis



Consumables None



Full DR











Maintenance Cost (\$)(\$)





CR



Film

Screen

CR



IP Barcoding



Diagnosis



40 sec.

60 sec.







RoseM's Benefits



Increase Image Quality



Increase no. of patients



Improvement of Mammo. PACS
Usability



Increased effectiveness of workflow



Increased space usability



Multiple applications [Mobile Examination]







Savings in CR maintenance cost



Reduction in Film/IP
Cassette maintenance and
human resource



RoseM Performance: Image Quality

76 µm the smallest pixel size



Micro-calcification detection is known to be a critical factor in the early detection of breast cancer. Visualization of micro-calcifications in a breast image is dependent on the performance and characteristics of the image receptor device (detector) such as the small pixel size.

Low noise and high DQE RoseM with its 76 µm holds the title for having **the smallest pixel size** amongst CsI (Indirect) mammography detectors. Also, with its low noise electronics, high DQE and MTF is achieved resulting in a sharp and highly defined mammography imaging optimized for screening and diagnostic purpose.


RoseM promises

- ✓ High quality images paralleling high-end FFDM
- ✓ Faster image acquisition
- ✓ Easier and more accurate diagnosis






RoseM Performance: Image Quality




RoseM Performance: Image Comparison

RoseM


vs. Full Digital Mammo System



RoseM



High End FFDM




Low- Mid End FFDM




RoseM Performance: Image Comparison (ACR Phantom)

RoseM


vs. Full Digital Mammo System




System A (a-Se)



System B (CsI)



RoseM







RoseM Performance: Image Comparison

RoseM

vs. CR Needle IP









RoseM Performance: Image Comparison

RoseM

vs. CR Needle IP









RoseM Core Technologies: TRUVIEWTM ART

Advanced image Reconstruction Technology




Image sharpness of an object in a conventional image is reduced due to light scattering. TRUVIEW® ART, DRTECH's unique reverse filtering technology reconstructs and improves image sharpness to increase the possibility detecting abnormalities including microcalcifications.






CONVENTIONAL IMAGE





RoseM Core Technologies: Lossless AED

RoseM's innovative Lossless AED function guarantees **faster and more reliable AED operation**. Conventional AED must have three steps of X-ray sensing, panel reset and charge accumulation. Integration time is a loss due to the extra period required to capture the x-ray sensing from immediately after x-ray shooting to panel reset.



When separate sensor modules within the detector system are used, x-ray loss is inevitable. With thick objects the loss rate can increase more, as panel reset and signal integration only occurs after the x-ray is detected. Lossless AED innovatively improved the reliability of sensitivity through optimization of panel driven scheme.



RoseM Core Technologies: RConsole1

Imaging Software for Mammography





Light dispersion is removed and image shapness increased to improve the detection of abnormalities including microcalcifications



9 image style options are provided to cater for various user preferences



With one click, image layout is easily adjusted



Patient's current and past image can be easily compared with instant image recall



Fast inspection of image quality is possible with large thumbnail display



Supported in 12 languages for convenient use worldwide




RoseM Proven Market Performance in Korea

RoseM has created extra demand for digital mammography in Korea!

Mammography Installation Status in Korea

* Total Annual Sales of New mammography: ~100units



Source: Korean Institute for Accreditation of Medical Imaging

Market Competition Dynamics

Unit: \$ 1,000

• Low: \$15~25

• Mid: \$35~50

• High: \$50~100

• Low: \$60~80

• Mid: \$90~150

• High: \$180+




RoseM is creating **NEW MARKET** while new system market size is maintained!




RoseM Proven Market Performance in Korea

RSM Installation by system manufacturer



RSM Installation by system model






RoseM Proven Market Performance in Korea

- The initial sales from **health screening hospitals** → gradually expanded to **breast surgery clinics**
- RSM is also being used in **specialized image acquisition** (breast magnification, breast needle localization, specimen, etc.)

RSM installation Sites in Korea






RoseM Marketing in Korean Market

Total Annual Market Size: ~ 100 units

Total RoseM Sales in 2016: **85 units**

How could this happen?



- Excellent Image Quality
- Easy Installation
- Versatile usages
 - For an old analog system
 - With a new analog system
 - For a mobile health bus

- Premium positioning (but still price competitive)
- Let Image Quality Speak
- Support with a value comparison chart (RoseM vs. CR / FFDM)

Aggressive product demonstrations

Right

- Get endorsements from Luminaries
- Effective Push & Pull

"We are comparable with FFDM, Not with CR!"



RoseM Marketing in Korean Market

Effective Push & Pull Marketing

PUSH

with consistent marketing messages

make 'Digital' something critical

PULL

Aggressive Demonstrations







Publicize & arouse the benefits of Digital Mammography

Direct Mails (Quarterly)









Remove customer uncertainty with publicity & awards

Conference & Seminars









Viral marketing with Luminaries



RoseM Marketing in Korean Market

Aggressive sales strategy - Demo. and Survey

Demo / Sales			//24300	Installation	on		
Information							
Type	Demo Installati	01	635	o Sales Installa	ton		
Requested Installation	2017, 551,00	20 11 MM, 00. X Destion		of the demonstration is 5 working days.			
Requested Uninstallation	20 YY, MM, CO. X Uninstal		distion can be discussed in special circumstances				
Sité Information							
Name	p 500+	0.400-4	90	o 300-399	p 200-299		
Site Bed Number	D 100-199	0.100+		Other (1 0 200 200		
Address		-			1		
Contact 1			Contact	Hospital	8		
Contact 2				Mobile			
System Information							
Manufacturer	is Hologic	n GE		o SEMENS	o MS Goto		
Model	o Planned	o Toshit	M.	o Metaltronica	DT3 tt		
Acadoton	to Full-Auto		o Semi-Auto		:: Manual		
Cassette Bucky Type	o MIN R2 type (f			□ MN S type ()			
Environment Information							
Imaging Tool	n CR scanner us			:: Analog film as	ė		
CR Scanner Brand	to Others (o Konic	Minota CR	G Carestream C	R GAGFACR		
LAN Port (mitro)	:: Present			© Not precent	1		
FACS Installation status	is PACS present	è		D PACS not pre	oent		
PACS Viewer Status	o Required			○ Not Required			
DICOM Penter Access Mammography Room	in Access	dended from syste	-	 Not present Separate shie 	14		
Setting	p Block-out area		No.	independent l			
ACR Phanton	a Present	7.5.		⊕ Not precent	70.000,000		
Diagnostic Monitor Monitor Desk	□ Present □ Required			© Not present			
SIONSOF LIPER	Dividence			o Not Required			
Online Support	100						
Support Requerment	to Year			o No			
Reason							
Extra Installation				rz No			
Examination Bus	p Yes						
	p Yes o Required			o Not Required			

	만경	도 설문	지					
	- 186	-		211				
		and the same	[원/진료과] 성 함	205-				
reconant v			9 0					
제품 정보	CONTRACTOR OF THE PARTY OF THE				_	_		_
당사 제품 정보대 대학 이란 때								
□ DM자료 □ 선은광교기사	D Set	[영구기사 [지원속생	□ 71€	8.0	_)	
. 데모·설치 만족도								
당사 매모 지원함의 설치에 안해	하십니까?				100	0	1 04	۵.
느 당사 매모 지원팀의 설치에	만족하지 않	요시다면 이유:	무엇입니	The r				
4						_		1
에므 기간에 만족하십니까?					0		04	₽ □
느 때로 기간에 만족하지 않으신	나면 적절			_	-			_
THE RESIDENCE OF THE PARTY OF T	3 14 S	⊕ 20 ¥	(5) 10)	발 또는 :	그 이상		ngersenne	
에모에 안족하십니까?					de	0	아니	<u> </u>
느 애모에 만족하지 않으시다면	U.S. V U.S.					-		
① 설치시간 ① 영상품질	① 사용성	④ 서비스대	8 57	타 (1		
IL 영상품질 만족도								
당사 제품의 영상 용장에 만족하	사십니까?				-78	0	94	20
나 영상 품질에 만족 또는 불만	축 하셨다면	이유는 무엇입	니까?				Total Co.	
① 명상 원목적	申申4	김 검을 정확성	0	디지털 양	상자	제 단 ⁴	투도:	
® 영상 파라메터 최착화	© 711 (1						
당사 제품의 영상 품질을 이게	각 성부사항	를 기준으로 5	병 만점으로	평가하	주세요			
 지방조직과 유선조직의 식별성 	4			- 1	2	3	4	5
- 지밀한 유선조직간의 식별능력	4			1	2	3	4	5
- 미세석회화 식별능력				1	2	3	4	5
- 종교의 식별능력				-1	2	3	4	5
- Artifact 발생 정도				11	2	3	14	5
RoseM 1824C 영상에 대한 의견	(교통) 보면		유롭게 기업	하우리의			444	
						_		
						_		
v. 구매의향			(ouzir		78	-	lau	0 11
V. 구매의향 . 당사 제품의 영상품질 확인 후.	제품 신뢰의	네 변화가 있었	-		_	_	-	_
V. 구매의향 - 당사 제품의 영상품질 확인 후, - 당사 제품의 영상품질 확인	제를 신뢰되 후 신키도아	대 변하가 있었 I 변화가 없거니	안족하지	않으셨다	_	_	-	_
V. 구매의향 당사 제품의 영상품질 확인 후, 노당사 제품의 영상품질 확인 - 다당사 제품의 영상품질 확인 - 다당사 제품의 명상품질 확인	제를 신화되 후 신체도에 용성 급 :	에 변화가 있었 [변화가 없거니 기준 장이의 노	한쪽하지 후화 ㅁ 7	일으셨다 다 (면 이4	\$ ± 4)	옛입니	項?
V. 구매의향 - 당사 제품의 영상품질 확인 후, - 도자 제품의 영상품질 확인 - 다 당사 제품의 영상품질 확인 - 다 영상품질 - 다 가격 - 다 사	제를 신뢰되 후 신뢰로야 용성 급 대했을 때 구	대 변하가 있었 I 면화가 없거니 기준 장미의 노크	안족하지 위화 급 / 있으십니까	일으셨다 다 (면 이4	_	옛입니	_



RoseM Testimonials

Testimonials by Korea's leading Breast Cancer Specialists and Radiologist



Dr. Yoo Hyun, Ahn

"We benefited greatly by using RoseM during our peak screening season. Our patients were much more happier and the screening time is shortened and the waiting time decreased. I am also very happy with the image quality. RoseM has helped greatly in increasing my confidence in making accurate diagnostic judgments"

Prior to converting to RoseM, we were using the CR system. When considering the conversion to digital mammography, we considered products of globally recognized system makers but decided to purchase RoseM as the image quality was very impressive and could be compared to that of more expensive full DR systems. My colleagues and I are very please with RoseM so far and we are especially happy that we are able to get high quality images at a very economical price"

Wonju Yeonse Min Clinic Radiologist Dr. Dong Ho Yeom

"I have a lot of patient that come for a second opinion after finding abnormalities at their initial screening. Therefore I place high emphasis on image quality. At first I had my doubts about an upgrade product but after using RoseM for the past 4 months, I am very pleased with the quality of images and I am recommending RoseM 1824C onto my colleagues"

Dr. Hyun Joo, Choi
Vice Chairman,
Korean Society of Breast Imaging



RoseM Q&A



Q. Is perfect digitalization possible simply with a RoseM detector?

Yes! Perfect digital upgrade is possible simply by inserting RoseM to your existing analog mammography device. RoseM is an officially approved medical device for mammography application with MFDS (Korea), CE (European) and FDA (USA) approval.

Q. Is the image quality of RoseM comparable with a FFDM?

Yes! RoseM comes along with RConsole1 image processing software designed by radiologists specifically for RoseM. RConsole1 is equipped with TRUVIEWTM ART, an image processing algorithm optimized for mammography imaging, which enhances the image sharpness. With this advanced technology, image quality paralleling high-end FFDM is possible and it has been proven with entrustment of many renowned breast cancer authorities in Korea choosing RoseM for their diagnostic practice.


Q. Is RoseM compatible with all analog mammography devices?

Yes! RoseM is compatible with all analog mammography units in operation worldwide. The performance of your mammography system will not be affected in any way by the addition of RoseM. The most important component of a mammography device is the detector. The performance of the detector mostly determines the overall image quality.

Q. Is it difficult to transit to RoseM from film-screen or CR?

No! Screening using RoseM is much easier and hassle free than using a film-screen or CR. By digitalizing with RoseM, you can eliminate one process (film printing and/or barcoding). In addition, full training is provided for users and our RConsole1 software will also help you further with its easy and efficient workflow.

RoseM Specifications



Conversion Layer

Pixel Pitch

Resolution

MTF

Preview Time

Cycle Time

Data Interface

AEC Sensing

Chest Wall Distance

Weight

Overall Dimension

Image Processing

Direct Deposition Csl

76 um

2,304 x 3,072

65 @ 3 lp/mm

≤4sec.

10 - 13 sec.

Giga Ethernet

System Compatible AEC

 \leq 2.0 mm

0.92 kg (Detector Only)

194.5 x 267.5 x 14.2 mm

TRUVIEW@ ART

Direct Deposition Csl

76 um

3,072 x 3,840

65 @ 3 lp/mm

≤4sec.

10 - 13 sec.

Giga Ethernet

System Compatible AEC

≤ 20 mm

1.2 kg (Detector Only)

253.7 x 327.5 x 14.2 mm

TRUVIEW@ ART



RoseM AEC Compatibility

Company	Brand name		AEC compatibility					
Company	Dianu name	Enable	Disable	Unknown				
Siemens	MAMMOMAT 1000	0						
	MAMMOMAT 3000 Nova	0						
	MAMMOMAT Select	0						
	Instrumentarium Performa	0						
	Alpha RT	0						
	Alpha ST	0						
GE	Senographe 600T	0						
	Senographe 700T	0						
	Senographe 800T	0						
	DMR Plus	0						
	Affinity	0						
Hologic	MIV	0						
(Lorad)	M IV Platinum	0						
	HAWK-2M	0						
Philips	MammoDiagnost	0						
Planmed	Nuance Classic	0						
	Sophie Classic	0						
Toshiba	MGU-1000A	0						
Giotto	Analog Giotto Image	0						
Metaltronica	Lilyum	0						
Genoray	MX300/500	0						
	MX600	0						
Bemems	Pinkview-AT/BT/RT	0						
Demeins	GENDEX MAMEX	0						
Bennett	Bennett Contour	0						
Definett	MXR	0						
DK	ELMA-T3	\circ		0				

Thank you



