

yuwell



Nebulizer

405A Air-compressing Nebulizer

Technical data:

Input: AC100 -220V, 50/60Hz 350mA

Output: DC12V, 1000mA

Max nebulizing rate: ≥ 0.2 ml/min

MMAD: $2.62 \mu\text{m}$ (Compliance with European Standard EN13544-1)

Max pressure: ≥ 0.1 MPa

Max air flow: ≥ 3.5 L /min

Operation noise: approx 57dB(A)

Operation time: 20minutes on, 40 minutes off

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Sales Center Address: Huanyuan East Road No.1, Xuzhuang Software Park, Nanjing, Jiangsu Province, P.R. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn

www.yuwell.com

f | t | in | i | Yuwell Global