

LITE MEDICS
Prime



to be **prime** is an attitude

www.litemedics.com - info.litemedics.com



OPEN PARAMETERS for expert users **who want the best**

no more cable with the **WIRELESS FOOTSWITCH**

EXTRA PROTOCOLS available to upgrade your practice day by day


Specific handpiece for **THERAPY, WHITENING** and **SURGERY**

up to **12W**



WAVELENGTH	980 NM
MODEL NAME	LITEMEDICS PRIME
LASER PEAK POWER	12 W
POWER SUPPLY INPUT	100 - 240 AC 50/60 HZ
PULSE LENGTH	SNP - CONTINUOUS
PULSE FEATURES	SUPER SHORT PULSE
FREQUENCY	CONTINUOUS - 25KHZ
DIMENSIONS	250X230X100 MM LXPXH
NET WEIGHT	2.5 KG
MEDICAL CLASS	II B
LASER CLASS	4

LITEMEDICS is a line of product made by LAMBDA Spa - via dell'Impresa 1, 36040 Brendola (VI) ITALY



CE
0476

AM ST 148