

# yuwell



MEDVAL  
healthcare



## Positive Airway Pressure Units

### YH450

#### Specification

Mode:CPAP、Auto CPAP

Pressure:4-20cmH<sub>2</sub>O

Efficient humidifier、Smart Start/Stop、FPS-Tech、

Leak Reminder、Altitude Adjustment、Leak Compensation、

Data Management、Central Apnea Detection、Smart Humidification

---

#### JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Sales Center Address: Huanyuan East Road No.1, Xuzhuang Software Park, Nanjing,  
Jiangsu Province, P.R. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn  
www.yuwell.com

f | t | in | i | Yuwell Global