

# yuwell



MEDVAL  
healthcare



## Positive Airway Pressure Units

### YH825

#### Specification

Mode:CPAP, S, T, ST, VGPS

Pressure:4-25cmH<sub>2</sub>O

Efficient humidifier, High Pressure Reminder, Leak Reminder,  
Asphyxiation Reminder, Trigger Adjustment, Cycle Adjustment,  
Slope Adjustment, Tidal Volume Adjustment, Altitude Adjustment,  
Leak Compensation, Data Management, Central Apnea Detection,  
Smart Humidification

---

#### JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Sales Center Address: Huanyuan East Road No.1, Xuzhuang Software Park, Nanjing,  
Jiangsu Province, P.R. China

Tel: 86-25-8713 6530 Fax: 86-25-8713 6530

E-mail: yuwellsales@yuyue.com.cn ; sales@yuyue.com.cn  
www.yuwell.com

f | t | in | @ | Yuwell Global