





5kW FPD Surgical C-arm **Garion**





Experience The Difference

Fine-Tuned Movement 5kW FPD Surgical C-arm



Garion is a compact and versatile digital surgery system. It is aimed to help you save space and improve work efficiency. It has minimal post-sales expenditures.

Garion for complete satisfaction

Garion is equipped with a high resolution FPD.

Its composition, along with other components such as positive/negative inversion, realtime zoom & rotation, as well as advanced noise reduction and optimized contrast/brightness, guarantees much clearer and enhanced image quality.

Attachable Grid

The Grid can be attached and detached, allowing low-dose examination for each part.

Dual Laser Pointer

Using the dual Laser pointer, it is possible to set the patient's position without unnecessary exposure.



Software

Powerful Digital Image Processing and Display

Whatever the surgical and clinical environment is, you will benefit from the most optimized output.

The combination of Dynamic FPD means that Garion provides the highest quality on screen image.

Our advanced image-processing technology enables you to utilize various functions such as image flip (horizontal & vertical), image rotation (360°), realtime zoom, noise reduction, and realtime edge enhancement, etc.













2. Specification

X-RAY GENERATOR

1) Power output	5.0 kW
-----------------	--------

2) Fluoroscopy

-Continuous Mode kV / mA range $40 \sim 110 kV / 0.3 \sim 6.3 mA$

-Pulsed mode kV / mA range 40 ~ 125 kV / 0.3 ~ 32 mA Frame rate 1, 2, 5, 10, 15 fps

3) Radiography kV/mAs range $40 \sim 125 kV / 0.1 \sim 100 mAs$

X-RAY TUBE

1) Anode Type Rotating Anode
2) Anode heat capacity 200kHU

3) Focal Spot 0.3 / 0.6 mm

DETECTOR

1) Sensor Panel Technology a-SI:TFT with direct dep. Csl

2) Active Area 209.9 x 209.9 mm 3) Matrix Size 1,024 x 1,024 pixels

4) Pixel Pitch 205 µm

GEOMETRY

1) Rotation Orbital 150°

Pivot 360° (±180°) Swing 30° (±15°)

Swing 30° (±15°)
2) Movement Vertical 200 mm

Horizontal 400 mm

3) Depth 640 mm 4) Free Space 784 mm

4) Free Space 784 mm 5) SID 980 mm

WORKSTATION

1) CPU Intel Core i5 6GEN

2) Memory 16 GB 3) Hard Disk 1 TB

4) UPS Power Supply Time 60 Min. (Load 90W)

MONITORS

1) Workstation Display 19" x 2

2) Operation Display 19" Touch

XX Note: Specifications subject to change without notice.



3. Service

Cutomer Oriented

Your values are of utmost importance to us. As a leading medical diagnostic device manufacturer, we take full responsibility in satisfying your needs by providing high quality products to the global market. Also, by collaborationg user friendly interface and design, we actualize your desires and needs. We promise to stand by your side and always strive to add pleasure at your work place.

Approved Qulaity

We put great deal of emphasis on well-being of mankind. SG Healthcare upholds this value by manufacturing reliable and reasonable medical diagonostic devices. And also by distributing our quality products all over the world through fortified SG Healthcare's global network, we wish to present better life for every little breath.

24/7 Fast Response

We embrace all challenges and counter act with pride and quickness to all claims. Fast response offers customers with satisfaction and improves reliavility of SG Healthcare. We strive to provide the most comprehensive service for every claim, regardless of location and time. This is how we have maintained friendship over 60 countries worldwide.



