

LAFOMED



8L | 12L



18L | 23L

Program	121 C solid	121 C Hollow	121 C Porous	134 C solid	134 C Hollow	134 C Porous	Prion	User-Defined	BD/Helixa Test	Vacuum Test	Clean
Vacuum Times	1	3	3	1	3	3	3	User-defined	3	1	3
Temperature (C)	121	121	121	134	134	134	134	User-defined	134	/	105
Pressure (KPA)	110	110	110	210	210	210	210	User-defined	210	-30	20
Time (Min.)	20	20	20	4	4	4	19	User-defined	3.5	15	5

Premium 8L	Premium 12L	Premium 18L	Premium 23L
Capacity	8L	Capacity	12L
Chamber Size	Φ17*32cm	Chamber Size	Φ20*36cm
Display	LCD	Display	LCD
Lock	Double lock	Lock	Double lock
Product Size	63*43*36cm	Product Size	63*43*36cm
G.W./N.W	52/46KGS	G.W./N.W.	55/47KGS
Power	220V/50Hz 1.95KW	Power	220V/50Hz 1.95KW
Output	USB+Printer	Output	USB+Printer
Capacity		Capacity	18L
Chamber Size		Chamber Size	Φ25*35cm
Display		Display	LCD
Lock		Lock	Double lock
Product Size		Product Size	62*50*43cm
G.W./N.W		G.W./N.W.	61/53KGS
Power		Power	220V/50Hz 1.95KW
Output		Output	USB+Printer
Capacity		Capacity	23L
Chamber Size		Chamber Size	Φ25*45cm
Display		Display	LCD
Lock		Lock	Double lock
Product Size		Product Size	66*50*43cm
G.W./N.W		G.W./N.W.	63/55KGS
Power		Power	220V/50Hz 1.95KW
Output		Output	USB+Printer

LAFOMED



Ningbo Haishu Life Medical Technology Co.,Ltd.

No.1 Jinghui Rd., Hengjie Town, Haishu, Ningbo, Zhejiang.
315181, China PRC

Tel: 0086 574 8700 6890 Fax: 0086 574 8828 3099

E-mail: sales@lafomed.com



LAFOMED



AUTOCLAVE STEAM STERILIZER

PREMIUM LINE

Prevent Infection. Create Safety.

Premium Line



B Type

In compliance with international standards



Built-in Printer And USB Port

Easy & safe data saving



**7 Preset Program, 2 Test Program
1 Custom Program**

Suitable for all kinds of instruments



Drying Class Adjustable

To achieve perfect drying



Water Quality Checking Module

Prevent steam generator and valves from getting blockage



Reliable Locking System

- Double lock (mechanical and electrical)
- Two position mechanical lock design



**Automatic Filling Option
With Ro Water System**



**Bestseller In
European Market**

