



MEDVAL
healthcare

star


RF Electrode



star RF Electrode

star RF Electrode-Fixed

- Internally Cooled electrode designed with a small lightweight ergonomic handle
- Wide selection of Needle Gauge, Needle Length and Active Tip
- Positive user feedback on sharper Needle Tip for easier percutaneous approach
- Cost effective and suitable for RFA in most organs



Gauge	Length	Active tip
19G	7cm	4mm, 5mm, 7mm, 10mm, 15mm, 20mm
	10cm	5mm, 7mm, 10mm, 15mm, 20mm
18G	7cm	4mm, 5mm, 7mm, 10mm, 15mm
	10cm	5mm, 7mm, 10mm, 15mm, 20mm
	15cm	5mm, 7mm, 10mm, 15mm, 20mm, 30mm
17G	7cm	4mm
	15cm	7mm, 10mm, 15mm, 20mm, 25mm, 30mm, 40mm
	20cm	10mm, 15mm, 20mm, 25mm, 30mm, 40mm
	25cm, 35cm	10mm, 20mm, 30mm, 40mm
15G	15cm	15mm, 20mm, 25mm, 30mm, 40mm
	20cm	20mm, 25mm, 30mm, 40mm
	25cm, 35cm	20mm, 30mm, 40mm

