

Backrest



Legrest



Auto Contour



Knee Break



- Electrical operated backrest and legrest adjustment
- Fowler and Vascular positions
- Hand remote control
- Lockable and detachable PP head and foot boards
- Auto-contour
- Easily cleanable, lockable and fold-away PP side rails
- Auto-regression system for backrest
- Easy removable ABS mattress platform
- Angle gauge for back.
- Electrostatic powder coated metal frame
- Footrest adjustment by ratchet
- Height adjustable stainless steel IV pole
- 125mm diagonal lockable castors
- Protective plastic corner bumpers
- Dual side manual CPR levers at backrest (optional)
- Trendelenburg position (optional)
- IPX6 standard electronic system

www.optium.com.tr

TECHNICAL DETAILS

- Overall Length : 216 cm
- Overall Width : 102 cm
- Height : 50 cm
- Backrest Angle : 0°-70°
- Legrest Angle : 0°-40°
- Castor Diameter : 12,5 cm
- Underbed Clearance : 12,5 cm
- Safe Working Load : 250 kg

ELECTRICAL DATA

- Power In : 1.5A max. At 230V AC 50-60 Hz
- Electrical Shock Protection : Class I, Type B
- Liquid Ingress Protection : IPX6



Accessories



X-Ray Translucent HPL Backrest




X-Ray Cassette Holder at Backrest



Rechargeable Battery Backup



Dual Sided Manual CPR at Backrest



Lifting Pole



Trendelenburg by gas spring



Bed Extension



Orthopedic Traction System

Color Options



P01



P02



P03



P04



P05



P06



P07



P08

www.optium.com.tr

