



Backrest



Legrest



Auto Contour



Knee Break



High / Low



- Electrical operated backrest, height and legrest adjustment
- Fowler and Vascular positions
- Hand remote control
- Lockable and detachable PP head and foot boards
- Auto-contour
- Easily cleanable, lockable and fold-away PP side rails
- Auto-regression system for backrest
- Easy removable ABS mattress platform
- Angle gauge for back
- Electrostatic powder coated metal frame
- Footrest adjustment by ratchet
- Height adjustable stainless steel IV pole
- Protective plastic corner bumpers
- 125mm diagonal lockable castors
- 125mm central lockable castors (optional)
- Dual side manual CPR levers at backrest (optional)
- Trendelenburg position by gas spring (optional)
- Bed Extension (optional)
- Battery backup (optional)
- IPX6 standard electronic system

### TECHNICAL DETAILS

• Overall Length	: 216 cm
• Overall Width	: 102 cm
• Height Range	: 37-80 cm
• Backrest Angle	: 0°-70°
• Legrest Angle	: 0°-40°
• Castor Diameter	: 12,5 cm
• Underbed Clearance	: 12,5 cm
• Product Weight	: 100 kg
• Safe Working Load	: 250 kg

### ELECTRICAL DATA

• Power In	: 1.5A max. At 230V AC 50-60 Hz
• Electrical Shock Protection	: Class I, Type B
• Liquid Ingress Protection	: IPX6



### Accessories



X-Ray Cassette Holder at Backrest



X-Ray Translucent HPL Backrest



Rechargeable Battery Backup



Central Brake Castors 150mm



Dual Sided Manual CPR at Backrest



Lifting Pole




Oxygen Cylinder Holder



Trendelenburg by gas spring



Bed Extension



Orthopedic Traction System

### Color Options



P01



P02



P03



P04



P05



P06



P07



P08

www.optium.com.tr

