



www.optium.com.tr



- Hydraulic height adjustment by dual sided pedals
- Manual backrest adjustment by gas spring
- Trendelenburg & Reverse Trendelenburg by gas spring
- Synthetic leather coated mattress
- Collapsible side rails
- Undercarriage storage for equipment
- Protective bumpers
- Swivel castors, 2 with brake

TECHNICAL DETAILS	
Width:	78 cm
Length:	215 cm
Trendelenburg Angle:	12°
Reverse Trendelenburg Angle:	12°
SleepSurface Dimensions:	60x190 cm
Height Range:	53- 87 cm
Castor diameter:	150cm

**OPTIUM MEDİKAL LTD. STI.**

Airport Industrial Zone Hidayet Street No: 23 Oğuzeli, GAZIANTEP, TURKEY  
Tel/Fax: + 90 850 807 62 84 e-mail: info@optium.com.tr web: www.optium.com.tr

